

# AD-BAND<sup>TM</sup>

**"BLOW THE BAND"** to achieve complete Hemostasis



Radial Artery Compression Device

### Designed to safely and comfortably stop bleeding

- Specially designed for radial artery hemostasis post percutaneous puncture.
- Transparency material ensures clear visual control of hemostasis procedure.
- Unilateral valve connected to balloon permits accurate pressure adjustment, providing personalized hemostasis strategy.
- Adjustable band length provides precise fit for patient comfort.

### Ordering information:

Catalogue No.	Balloon structure	Bandage length (cm)	Max injection air volume (ml)
ADB 24	2	24.5	20



Changing lives for the better...

### ADVANCED LIFE SCIENCES PVT. LTD.

D-22, Okhla Industrial Area, Phase-I, New Delhi-110020, INDIA

Tel. : +91 (0)11 - 26812108, 26371127, 26812107

Fax : +91 (0)11 - 26371128

Email: [customerservice@alspl.com](mailto:customerservice@alspl.com); Website: [www.alspl.com](http://www.alspl.com)