


Ad-DucerTM


Introducer Sheaths

Color Code
4F
5F
6F
7F
8F
9F
10F



Femoral Sheaths with Dialator

Reference Number	Sheath Size	Sheath Length (cm)	Guidewire Size (J Tip)	Introducer Needle
AFIS 4F-9-021IN	4F	9	0.021"x45cm	20Gx45cm
AFIS 5F-11/23-038IN	5F	11/ 23	0.038"x45cm	18Gx7cm
AFIS 6F-11/23-038IN	6F	11/ 23	0.038"x45cm	18Gx7cm
AFIS 7F-11/23-038IN	7F	11/ 23	0.038"x45cm	18Gx7cm
AFIS 8F-11/23-038IN	8F	11/ 23	0.038"x45cm	18Gx7cm
AFIS 9F-11/23-038IN	9F	11/ 23	0.038"x45cm	18Gx7cm
AFIS 10F-11/23-038IN	10F	11/ 23	0.038"x45cm	18Gx7cm



ADVANCED LIFE SCIENCES PVT. LTD.

D-22, Okhla Industrial Area, Phase-I, New Delhi-110020, INDIA

Tel. : +91 (0)11 - 26812108, 26371127/28


Mob. : +91 9818237529, 8920716599

Email: customerservice@alspl.com; Website: www.alspl.com

Ad-Ducer™


Transradial Sheaths with Dialator

Reference Number	Sheath Size	Sheath Length (cm)	Hydrophilic wire size (Straight Tip)	I.V. Catheter
ATIS5F-07/011-025	5F	7/ 11	0.025"x45cm	20Gx4.5cm
ATIS6F-07/011-025	6F	7/ 11	0.025"x45cm	20Gx4.5cm
ATIS7F-07/011-025	7F	7/ 11	0.025"x45cm	20Gx4.5cm



Radial Sheaths with Dialator

Reference Number	Sheath Size	Sheath Length (cm)	Guidewire Size (J Tip)	I.V. Catheter
ARIS5F-07-025	5F	7	0.025"x45cm	20Gx4.5cm
ARIS6F-07-025	6F	7	0.025"x45cm	20Gx4.5cm
ARIS7F-07-025	7F	7	0.025"x45cm	20Gx4.5cm



Peelable Sheaths with Dialator

Reference Number	Sheath Size	Sheath Length (cm)	Hydrophilic wire size (Straight Tip)	Introducer Needle
ADP7F-14	7F	14	0.038"x60cm	18Gx7cm
ADP8F-14	8F	14	0.038"x60cm	18Gx7cm
ADP9F-14	9F	14	0.038"x60cm	18Gx7cm

***All kits are customizable on request**

ADVANCED LIFE SCIENCES PVT. LTD.

D-22, Okhla Industrial Area, Phase-I, New Delhi-110020, INDIA

Tel. : +91 (0)11 - 26812108, 26371127/28

Mob. : +91 9818237529, 8920716599

Email: customerservice@alspl.com; Website: www.alspl.com

