



AD-Valve™

HEMOSTASIS VALVE SET



For total control and clean procedure.

AD-VALVE Hemostasis valve allows single-handed manipulation to open and close the valve while introducing devices into the vasculature. The hemostasis seal minimizes blood loss while giving the clinician the ability to manipulate interventional devices.



Changing lives for the better...




Ordering information:

Catalog No	Description
AVR	9F Rotating Hemostasis Valve
AVR-20	9F Rotating Hemostasis Valve with extension line and stop cock
AVRK	9F Rotating Hemostasis Valve, Insertion tool, torque device
AVRK-20	9F Rotating Hemostasis Valve with extension line and stop cock, Insertion tool, torque device




Ordering information:

Catalog No	Description
AVP	7F Pull-Push Hemostasis Valve
AVP-20	7F Pull-Push Hemostasis Valve with extension line and stop cock
AVPK	7F Pull-Push Hemostasis Valve, Insertion tool, torque device
AVPK-20	7F Pull-Push Hemostasis Valve with extension line and stop cock, Insertion tool, torque device



Ordering information:

Catalog No	Description
AVC	7F Click Hemostasis Valve
AVC-20	7F Click Hemostasis Valve with extension line and stop cock
AVCK	7F Click Hemostasis Valve, Insertion tool, torque device
AVCK-20	7F Click Hemostasis Valve with extension line and stop cock, Insertion tool, torque device



Ordering information:

Catalog No	Description
AVDR	9F Double Y Rotating Hemostasis Valve
AVDRK	9F Double Y Rotating Hemostasis Valve, Insertion tool, torque device



Ordering information:

Catalog No	Description
AVIT	20G Insertion Tool 100mm



Ordering information:

Catalog No	Description
AVTD	Torque Device 0.010"-0.016"



Ordering information:

Catalog No	Description
GTDXL	Glider torque device

Not for sale in USA



Changing lives for the better...

ADVANCED LIFE SCIENCES PVT. LTD.

D-22, Okhla Industrial Area, Phase-I, New Delhi-110020, INDIA

Tel. : +91 (0)11 - 26812108, 26371127, 26812107

Fax : +91 (0)11 - 26371128

Email: customerservice@alspl.com; Website: www.alspl.com