



## ARROW PERCUTANEOUS SHEATH INTRODUCERS

Another effective choice

### AT A GLANCE

- highly kink-resistant Polyurethane shaft
- reduced kinking at hub through stress-relief hub sleeve
- reduced peel-back effect due to special design
- protection against blood contamination through Raulerson Syringe
- hydrophilic coating on distal sheath enhances placement

Through the ARROW brand, Teleflex provides an immense assortment of high-end products for vascular access. Thanks to the outstanding construction, the sheath introducers offered here reduce the potential for complications in both arterial and venous applications.

An outstanding kink-resistance due to use of our blend Polyurethane and our special stress-relief hub sleeve at the hub; reduced instances of peel-back or buckling thanks to the thicker tip and the heavy-walled design; ease of use and added protection against blood contamination through our exclusive Raulerson Syringe ensure a safe performance throughout the whole procedure.

Moreover, the wide array of sheaths sets and individual sheaths given, assures that you will have the flexibility in choices for multiple vascular applications.



\* U.S. Patent Nos. 4,813,938, 5,045,065, and 5,913,814; additional international patents.  
Advanced Patent Nos. 6,477,402 and 5,484,419.

## ORDERING INFORMATION

"CW" SETS PSI				ARROW
REF.	SIZE	LENGTH	COLOR CODED HUB	QTY
CW-08403*	4 Fr.	7.5 cm	● pink	10
CW-08503*	5 Fr.	7.5 cm	● grey	10
CW-08503-A	5 Fr.	11 cm	● grey	10
CW-08603	6 Fr.	11 cm	● green	10
CW-08703	7 Fr.	11 cm	● orange	10
CW-08803	8 Fr.	11 cm	● blue	10
CW-08903	9 Fr.	11 cm	● black	10

"CP" SETS PSI				ARROW
REF.	SIZE	LENGTH	COLOR CODED HUB	QTY
CP-08403*	4 Fr.	7.5 cm	● pink	10
CP-08503*	5 Fr.	7.5 cm	● grey	10
CP-08503-A	5 Fr.	11 cm	● grey	10
CP-08603-P*	6 Fr.	7.5 cm	● green	10
CP-08603	6 Fr.	11 cm	● green	10
CP-08703	7 Fr.	11 cm	● orange	10
CP-08803	8 Fr.	11 cm	● blue	10
CP-08903	9 Fr.	11 cm	● black	10

"CL" SETS PSI				ARROW
REF.	SIZE	LENGTH	COLOR CODED HUB	QTY
CL-08403*	4 Fr.	7.5 cm	● pink	10
CL-08403-A	4 Fr.	11 cm	● pink	10
CL-08503*	5 Fr.	7.5 cm	● grey	10
CL-08503-A	5 Fr.	11 cm	● grey	10
CL-08603	6 Fr.	11 cm	● green	10
CL-08703	7 Fr.	11 cm	● orange	10
CL-08803	8 Fr.	11 cm	● blue	10
CL-08903	9 Fr.	11 cm	● black	10

### EACH "CW" SET INCLUDES:

- 1 – ARROW Radiopaque Polyurethane Sheath with Integral Side Port/Hemostasis Valve and attached 3-Way Stopcock
- 1 – Arrow Raulerson Syringe
- 1 – Vessel Dilator with Snap-Lock Feature
- 1 – 17 3/4" (45 cm) long dual purpose spring-wire guide (straight soft tip on one end, "J" tip on the other end) with Arrow Advancer™. Pediatric set (\*) has .021" dia. and all others have .035" dia. wire guide
- 1 – 18 Ga. x 2 1/2" T.W. Needle

### EACH "CP" SET INCLUDES:

- 1 – ARROW Radiopaque Polyurethane Sheath with Integral Side Port/Hemostasis Valve and attached 3-Way Stopcock
- 1 – Vessel Dilator with Snap-Lock Feature
- 1 – 17 3/4" (45 cm) long dual purpose spring-wire guide (straight soft tip on one end, "J" tip on the other end). Pediatric set (\*) has .021" dia. and all others have .035" dia. wire guide

### EACH "CL" SET INCLUDES:

- 1 – ARROW Radiopaque Polyurethane Sheath with Integral Side Port/Hemostasis Valve and attached 3-Way Stopcock
- 1 – Vessel Dilator with Snap-Lock Feature

### STANDARD (11 CM) INTRODUCER SETS INCLUDE:

- Sheath with 11 cm (4 3/8") useable length
- Dilator which accepts .035" to .038" guide wires
- Sheath with hydrophilic coating on 1" of tip
- Dilator with 5cm exposed length when snap-locked into sheath

### SHORT (7.5 CM) INTRODUCER SETS INCLUDE:

- Sheath with 7.5 cm (2 15/16") useable length
- Dilator which accepts .021" to .025" guide wires
- Sheath with hydrophilic coating on 1" of tip
- Dilator with 3cm exposed length when snap-locked into sheath

Caution: U.S. Federal law limits this device to sale by or on order of a physician.

#### DISTRIBUTED BY:

TELEFLEX MEDICAL HEADQUARTER EMEA, IRELAND · IDA Business & Technology Park · Athlone · Co Westmeath

Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.int@teleflex.com

UNITED KINGDOM Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com

SOUTH AFRICA Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com

WWW.TELEFLEX.COM

94 10 44 - 00 00 01 · REV A · 08 12 02 · All data current at time of printing. Subject to technical changes without further notice.