



Our Mission

Certificates/Awards

Company History

Contact Info

D'VILL™ Introducer

Long-Length, Kink-Resistant Hemostasis Introducer



For Introduction of Large Diameter Interventional Devices
Compatible with CP Stent® and BIB® catheters

Available in 12F and 14F Sheath Sizes - 30cm, 65cm, and 85cm



- **Hemostasis Valve**
Designed to minimize blood loss
- **Side Port and 3-Way Stopcock**
For convenient flushing
- **Braided Sheath**
For optimal flexibility and kink resistance
- **Radiopaque Marker Band**
For easy imaging

Indications for Use: Recommended for introduction of balloons, catheters, and other diagnostic and interventional devices.

NEW sizes highlighted

D'VILL™ Specifications

Introducer Size (FR)	Usable Length (CM)	Guide Wire(Inches)	Catalog No.
12	30	0.035	DV1230
14	30	0.035	DV1430
12	65	0.035	DV1265
14	65	0.035	DV1465
12	85	0.035	DV1285
14	85	0.035	DV1485

**-Each D'VILL set includes (1) one Sheath and (1) one Dilator
-0.035" guidewire compatibility**

The devices listed have CE mark approval. Any customer or distributor outside of the EU will be subject to their individual country's regulations in regards to the importation and sale of these products.

This publication could include technical or other inaccuracies or typographical errors. Changes are periodically added to the information herein; these changes will be incorporated in the new editions of the publication. NuMED Incorporated may make improvements and/or changes in the product(s) and/or the program(s) described in this publication at any time.