

Jobst Farrowwrap LBL

JOBST® FarrowWrap® Strong

Made from durable fabric with a soft inner layer. Recommended for patients with **moderate to severe edema**.

- Suitable for minor skin folds and shape distortion
- Reliable hold for fluctuating edema
- Liner included



Colors



Beige

Styles



Footpiece^{1,2,3}



Legpiece^{1,2,3}



Kneepiece¹



Thighpiece^{1,2,*}

1: RTW (Ready-To-Wear) | **2: CM** (Custom-Made) | **3: TTF** (Trim-To-Fit)

* Thighpiece must be worn with Kneepiece

Compression Class

30-40 mmHg

Garment Care



Machine
wash 30°



Do not
bleach



Do not
tumble dry



Lay flat
to dry



Do not
iron



Do not
dry clean

Jobst Farrowwrap LBL

JOBST

Compression for
**Empowering
Self-Care**



JOBST® FarrowWrap®
Adjustable Compression Device

Jobst Farrowwrap LBL

JOBST® FarrowWrap®

Easy-to-attach VELCRO® brand tabs support easy donning and doffing.



Easy Donning

Short-stretch technology provides low resting and high working pressures to enhance lymphatic and venous return.



Short Stretch Technology

VELCRO® brand hook & loop system allows patients to respond better to fluctuating edema by simply tightening or loosening the wrap.



Adjustable

Discreet design that's less bulky than multi-layer bandaging and can be worn under normal clothing.



Low Profile Design

* Overlapping bands not on JOBST FarrowWrap 4000
² Observational Study JOBST® FarrowWrap® 2019, data on file