

JOBST 

Compression for
Everyday
Style and
Comfort



JOBST[®] Opaque
Medical Compression Stockings
for Day-to-Day Comfort



JOBST® Opaque

Discreetly opaque, clearly comfortable

Help your patients lighten their daily lives with JOBST Opaque. Intelligently designed with high-quality multi-fiber yarns, JOBST Opaque stockings combine outstanding functionality with an especially soft and pleasant feel. Their refined opacity discreetly hides blemishes while effective moisture management and ease of donning mean that your patients can enjoy excellent comfort **all day, every day.**



Refined opacity and fine knit structure **discreetly cover** small imperfections



Effective moisture wicking **helps to keep feet dry** for all-day comfort



A special blend of **soft yarns** ensure comfortable compression



Quality and durability with **reinforced heel and toe zones**, designed for long-lasting compression



High performance moisture management

Keeping feet dry is an essential part of ensuring comfort. This is especially important for stockings that are intended for all-day use.

That's why an advanced moisture management system is integrated into JOBST Opaque stockings that utilizes **functional fibers and high-performance yarns**. This proven mix of yarns actively **wicks moisture away from the body** and through the garment where it quickly evaporates.

SoftFit - compression that stays in place all day long

Bring wearing comfort to a new level! JOBST Opaque Knee High stockings are optionally available with the **advanced SoftFit technology**. For the SoftFit top band additional silicone yarns are knitted in around the knee welt that **grip softly to keep the stocking up** without increasing tension or digging into the skin. Recommend SoftFit top bands to help provide secure fit and a comfortable feel for your patients.



Compression Classes*

CCL 1: 18 - 21 mmHg

CCL 2: 23 - 32 mmHg

Colors*



Natural

Styles*




Knee High



Thigh High

Measurements Table

	I	II	III	IV	V	VI
cG standard	43 - 48 cm	45 - 52 cm	49 - 56 cm	53 - 60 cm	56 - 64 cm	60 - 68 cm
cG wide/weit	49 - 57 cm	53 - 62 cm	57 - 67 cm	61 - 72 cm	65 - 77 cm	69 - 81 cm
cG volant	43 - 57 cm	45 - 62 cm	49 - 67 cm	53 - 72 cm	56 - 77 cm	60 - 81 cm
cG	43 - 57 cm	45 - 62 cm	49 - 67 cm	53 - 72 cm	56 - 77 cm	60 - 81 cm
cD	27 - 33 cm	29 - 36 cm	32 - 39 cm	34 - 42 cm	36 - 45 cm	38 - 48 cm
cC	28 - 34 cm	30 - 37 cm	33 - 40 cm	35 - 43 cm	37 - 46 cm	39 - 49 cm
cB1	23 - 27 cm	24 - 29 cm	26 - 32 cm	29 - 35 cm	31 - 37 cm	33 - 39 cm
cB	18 - 20 cm	20 - 22 cm	22 - 24 cm	24 - 26 cm	26 - 28 cm	28 - 30 cm
Maxi option	cT > 80 cm / cH > 95 cm	cT > 90 cm / cH > 105 cm	cT > 110 cm / cH > 125 cm		cT > 120 cm / cH > 135 cm	



LK Petite: 65 - 75 cm / Regular: 75 - 85 cm | LG Petite: 60 - 69 cm / Regular: 70 - 80 cm | LD Petite: 34 - 39 cm / Regular: 39 - 45 cm

* Please refer to your local JOBST website and/or local customer service representative for product assortment availability in your region.