

# JOBST®

## RELIEF

Medical Compression  
Stockings



- Quality and efficacy at a moderate price
- Available in a variety of styles and compression levels

**No 1**  
PHYSICIAN  
RECOMMENDED  
in the U.S.



JOBST RELIEF IS A  
**MEDICAL COMPRESSION STOCKING**  
IN A UNISEX DESIGN DELIVERING GRADUATED COMPRESSION



15-20 mmHg\*

20-30 mmHg\*




Knee High




Thigh High

	S	M	L	XL
<b>cC</b>	<b>28 - 38 cm</b> (11"-15")	<b>30 - 42 cm</b> (11 <sup>7</sup> / <sub>8</sub> "-16 <sup>1</sup> / <sub>2</sub> ")	<b>32 - 46 cm</b> (12 <sup>1</sup> / <sub>2</sub> "-18 <sup>3</sup> / <sub>8</sub> ")	<b>34 - 50 cm</b> (13 <sup>3</sup> / <sub>8</sub> "-19 <sup>5</sup> / <sub>8</sub> ")
<b>cB</b>	<b>18 - 21 cm</b> (7"-8 <sup>1</sup> / <sub>4</sub> ")	<b>21 - 25 cm</b> (8 <sup>3</sup> / <sub>8</sub> "-9 <sup>7</sup> / <sub>8</sub> ")	<b>25 - 29 cm</b> (10"-11 <sup>3</sup> / <sub>8</sub> ")	<b>29 - 33 cm</b> (11 <sup>1</sup> / <sub>2</sub> "-13")



	S	M	L	XL
<b>cG</b>	<b>40 - 62 cm</b> (15 <sup>3</sup> / <sub>4</sub> "-24 <sup>3</sup> / <sub>8</sub> ")	<b>46 - 70 cm</b> (18 <sup>1</sup> / <sub>8</sub> "-27 <sup>1</sup> / <sub>2</sub> ")	<b>54 - 78 cm</b> (21 <sup>1</sup> / <sub>4</sub> "-30 <sup>3</sup> / <sub>4</sub> ")	<b>60 - 81 cm</b> (23 <sup>5</sup> / <sub>8</sub> "-32")
<b>cC</b>	<b>28 - 38 cm</b> (11"-15")	<b>30 - 42 cm</b> (11 <sup>7</sup> / <sub>8</sub> "-16 <sup>1</sup> / <sub>2</sub> ")	<b>32 - 46 cm</b> (12 <sup>1</sup> / <sub>2</sub> "-18 <sup>3</sup> / <sub>8</sub> ")	<b>34 - 50 cm</b> (13 <sup>3</sup> / <sub>8</sub> "-19 <sup>5</sup> / <sub>8</sub> ")
<b>cB</b>	<b>18 - 21 cm</b> (7"-8 <sup>1</sup> / <sub>4</sub> ")	<b>21 - 25 cm</b> (8 <sup>3</sup> / <sub>8</sub> "-9 <sup>7</sup> / <sub>8</sub> ")	<b>25 - 29 cm</b> (10"-11 <sup>3</sup> / <sub>8</sub> ")	<b>29 - 33 cm</b> (11 <sup>1</sup> / <sub>2</sub> "-13")



JOBST®,  
an Essity brand

06/2019

Essity India Private Limited.  
E-902, Lotus Corporate Park, Jay Coach Junction,  
Opp. Raheja Ridgewood, Goregaon (East), Mumbai - 400063.  
Tel. : +91-22-6155 4300 | email : infoindia@essity.com  
©Registered Trade Mark. License User of Trade Mark.

