



# Surgifold™

Manifolds for trusted control



Clear poly carbonate body with multiple styles  
and pressure ratings.

Leak free smooth rotation with low torque handles.



Changing lives for the better...

# Large bore inner lumen throughout, including rotator


## Ordering information:

**Description :** MANIFOLDS - Two Port 500 PSI / 250 PSI



Catalog No	Description
SF202N-R	*ON* Right-Handed w/Rotating Adapter, 2 Port, 250 PSI
SF202F-R	*OFF* Right-Handed w/Rotating Adapter, 2 Port, 250 PSI
SF502N-R	*ON* Right-Handed w/Rotating Adapter, 2 Port, 500 PSI
SF502F-R	*OFF* Right-Handed w/Rotating Adapter, 2 Port, 500 PSI

**Description :** MANIFOLDS - Three Port 500 PSI / 250 PSI



Catalog No	Description
SF203N-R	*ON* Right-Handed w/Rotating Adapter, 3 Port, 250 PSI
SF203F-R	*OFF* Right-Handed w/Rotating Adapter, 3 Port, 250 PSI
SF503N-R	*ON* Right-Handed w/Rotating Adapter, 3 Port, 500 PSI
SF503F-R	*OFF* Right-Handed w/Rotating Adapter, 3 Port, 500 PSI

**Description :** MANIFOLDS - Four Port 500 PSI / 250 PSI



Catalog No	Description
SF204N-R	*ON* Right-Handed w/Rotating Adapter, 4 Port, 250 PSI
SF204F-R	*OFF* Right-Handed w/Rotating Adapter, 4 Port, 250 PSI
SF504N-R	*ON* Right-Handed w/Rotating Adapter, 4 Port, 500 PSI
SF504F-R	*OFF* Right-Handed w/Rotating Adapter, 4 Port, 500 PSI

**Description :** High Pressure 3 way stop cock 500PSI / 1200 PSI



Catalog No	Description
SFHP3WFL 500	3 way stop cock 500PSI fixed luer lock
SFHP3WFL 1200	3 way stop cock 1200PSI fixed luer lock
SFHP3WRLL 500	3 way stop cock 500PSI rotating luer lock
SFHP3WRLL 1200	3 way stop cock 1200PSI rotating luer lock

Not for sale in USA



Changing lives for the better...

**ADVANCED LIFE SCIENCES PVT. LTD.**

D-22, Okhla Industrial Area, Phase-I, New Delhi-110020, INDIA

Tel. : +91 (0)11 - 26812108, 26371127, 26812107

Fax : +91 (0)11 - 26371128

Email: customerservice@alspl.com; Website: [www.alspl.com](http://www.alspl.com)